

EcoSac 5mm Tissue Retrieval System

from



Espiner Medical



In pediatric surgery, space is always at a premium.

JustRight Surgical is the exclusive distributor of the pediatric sized EcoSac 5mm Tissue Retrieval System from Espiner Medical.

EcoSac 5mm Tissue Retrieval System:

- Rip-stop nylon bag coated with polyurethane - prevents leaking
- EcoSac becomes translucent when wet for greater visibility
- Easily inserts through a standard 5mm cannula
- 60ml and 80ml bag sizes will fully deploy in the patient cavity
- Securely closes with a suture drawstring
- Easily manipulated and retrieved with a standard 5mm grasper

This specimen retrieval system comfortably fits down any 5mm cannula and provides durability, flexibility, and ease of use.



*This latest product offering from **JustRight Surgical** is another example of our commitment to providing “right sized” instrumentation for pediatric surgery.*



Product Name	EcoSac 5mm Tissue Retrieval System 60ml
Catalog Number	EMP60ECO-R-10
Package	Case of 10 retrieval bags
Shaft Length	45 cm
Shaft Diameter	5 mm
Volume	60 ml

Product Name	EcoSac 5mm Tissue Retrieval System 80ml
Catalog Number	EMP80ECO-R-10
Package	Case of 10 retrieval bags
Shaft Length	45 cm
Shaft Diameter	5 mm
Volume	80 ml

**THE FUTURE OF
SURGERY IS HERE,
AND IT'S TINY.**

