


JustRight™ 5mm Stapler



The Advantages of Thinking Small

At less than half the size of all others, Bolder Surgical's JustRight™ Stapler is the world's only 5mm stapling instrument, providing advantages for surgeons and patients.*

Access:


5mm shaft can get into tighter spaces than 12mm staplers.

Options for Port Placement:

5mm shaft provides flexibility in stapler placement allowing for optimal approaches to the surgical site.


Visualization:

5mm jaw takes up less of the visible field in the scope than does larger stapling devices allowing better visualization of the surgical field.



Provides Secure Staple Line

Classic B formation of 2.0mm titanium staples in 4 rows of 6 staggered staples provides secure staple line.




5mm Incision Site

Patient Benefits of 5mm incision vs. 12mm incision:

Patients may experience:

- Minimal post-op pain^{1,2}
- Minimal scarring²
- Minimal risk of port site herniation^{1,3,4}
- Faster recovery times¹
- Decreased potential for surgical site infection^{1,2}

JustRight™ 5mm Stapler & Reload



* Only stapler on the market with classic B-shaped staples that fits down any conventional 5mm trocar

References

1. Trocar Site Hernia. Hitoshi Tonouchi, MD, PhD; Yukinari Ohmori, MD; Minako Kobayashi, MD, PhD; Masato Kusunoki, MD, PhD, Arch Surg. 2004;139(11):1248-1256.
2. Incisions do not simply sum. Blinmin T. Surgical Endoscopy. 2010; 24(10):1746-52.
3. Trocar Site Hernia After Laparoscopic Cholecystectomy Uslu, H, Erkek, A, Cakmakm A, Kepenekci, I, Sozener, U, Kocaay, . Journal of Laparoscopic and Advanced Surgical Techniques, Volume 17, Number 5, 2007
4. Port-site incisional hernia: A case series of 54 patients. Lambertz, B.O. Stüben, B. Bock, R. Eickhoff, A. Kroh, C.D. Klink, U.P. Neumann, C.J. Krones . Annals of Medicine and Surgery 14 (2017) Be11

Ordering Information

JustRight™ 5mm Stapler	Catalogue No.	Staple line Length	Staple Leg Length	Shaft Length	Shaft Diameter
	JR-ST25-2.0-6	25 mm	2.0 mm	20 cm	5 mm
JustRight™ 5mm Reload	Catalogue No.	Staple line Length	Staple Leg Length	Shaft Length	Shaft Diameter
	JR-REL25-2.0-12	25 mm	2.0 mm	20 cm	5 mm

BOLDER
SURGICAL

331 S. 104th Street Suite 200, Louisville, CO 80027
P: 720-287-7133 F: 720-287-7135
www.boldersurgical.com

© July 2019, Bolder Surgical Holdings, Inc.